


## ECONOMICAL LATEX-FREE TOURNIQUETS



PMD's disposable tourniquets are made from high quality nitrile and contain no latex. We use a proprietary blend of materials that gives the performance of latex but is completely latex-free. Our color bands are pre-cut for immediate use. A great value!

Item Number	Description
970001	Blue non-latex tourniquet
970001/Y	Yellow non-latex tourniquet
970001/PR	Purple non-latex tourniquet
970001/G	Green non-latex tourniquet
970001/O	Orange non-latex tourniquet
970001/P	Pink non-latex tourniquet

Measurements: 1" wide x 18" long x .025" thick

Packing: 250 pieces/box, 1,000 pieces/case

**FOR MORE INFORMATION AND TO ORDER:**

**CALL US AT (732) 447-2587 OR EMAIL US AT SALES@PMDMFG.COM**